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IMPLA -

Implantatsystem Cylindrical

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Instructions

IMPLA -**Implant Systems Cylindrical**

IMPLA implants are made of pure titanium (ISO 5832-2). The IMPLA system comprises surgical, prosthetic and laboratory components and instruments.

Important notes, warning

IMPLA implants are part of a complete system and may only be used with the original components and instruments of the system. The use of components which are not part of this system voids the guarantee and rules out parts being replaced. The IMPLA system may only be used by dentists, doctors and dental technicians familiar with dental implantology including diagnosis, preoperative planning, surgical techniques and prosthodontics.

Consultation is done in written form or by electronic media as well as with demonstrations during special workshops. All information is based on the state of the art and latest scientific findings known at the time of launching the product. However, we urge the users of our products to check the suitability of the product in each case. This is especially important for techniques not expressly recommended by the manufacturer. In case of doubt, we encourage our customers to confer with Schütz Dental GmbH/Section IMPLA.

Reuse and sterilization

All abutments and components are for single use only and must not be reprocessed. Any reuse involves the risk of cross-contamination! All instruments, abutments and components are packaged non-sterile and must be sterilized (EN 13060) prior to placement (see page 26–35). Using unsterile components increases the infection risk for the patient.

Indications

IMPLA implants Cylindrical Implants are indicated for:

- restoring individual tooth gaps
- bridgework spanning interdental and free-end gaps, in accordance with conventional and special principles of prosthetic treatment
- restoration of edentulous mandibles and maxillae

Further indications:

Implants with a diameter of 3.6 mm

These implants are an alternative in case of limited jaw ridge width. They are, because of their lower mechanical strength when compared with the standard diameter (4.5 mm), intended for the following indications:

- a) restoring individual anterior teeth of the mandible
- b) restoring individual lateral anterior teeth of the maxilla
- c) four interforaminal implants to restore the edentulous mandible with a bar construction

Implants with a diameter of 4.5 mm

The IMPLA standard implants with a diameter of 4.5 mm are used for all indications of oral enossal implants for mandibles and maxillae. A prerequisite is an adequate minimal amount of bone base with a width of about 5.8 mm.

Implants with a diameter of 5.5 mm

These implants are indicated for the use in the posterior tooth areas of mandible and maxilla (especially for the restoration of molars). The jaw ridge width should not be less than 6.8 mm.

The implants are primarily indicated for delayed implantations.

The cylindrical implants with a length of 6.5 mm have been developed to offer an alternative to augmentative measures in cases with little vertical bone availability or atrophic bone.

Contraindications

When selecting patients, the general contraindications for dental surgery must be taken into account.

These include: Vascular diseases, uncontrolled diabetes mellitus and other metabolic or systemic diseases affecting bone regeneration, blood coagulation disorder, anticoagulative therapy, bone metabolism diseases.

Temporary contraindications: Acute and chronic inflammations, inadequate soft tissue present, chemotherapy or radiotherapy.

Oral contraindications:

- patient not really prepared to undergo total oral rehabilitation
- inadequate oral hygiene
- untreated parafunctions, e.g. bruxism
- inadequate amount of bone available
- inadequate interocclusal distance
- inadequate occlusion and/or articulation

Absolute contraindication

- Patients whose systemic condition does not allow for an implantation.

Important notes

When working with our implants and instruments intraorally, always take increased care to protect against aspiration! Primary stability directly after the insertion of the implant is absolutely necessary for a successful healing. During the healing phase, IMPLA implants are generally never exposed to a masticatory load.

Special information regarding implants with a length of 6.5 mm

In cases with a low residual crestal bone height, cylindrical implants with a length of 6.5 mm must only be used with great care and after a special assessment of the medical case.

In principle, cylindrical implants with a length of 6.5 mm are not suited to be used for individual tooth restorations or for supporting a closed prosthesis with no other support.

Surgical procedure

To create optimum conditions for the implant to heal, it is essential to treat the hard and soft tissues gently. It is advisable to use the countersink and bone thread cutter in order to ensure that the implant can be placed atraumatically applying as little torque as possible (maximum 30 Ncm). The so-called „Dense Bone Drills“ have been designed especially for the cortical bone preparation. They are meant to reduce the insertion torque for longer implants. As thermal trauma may impede or prevent the implant healing, heat generation should be kept to a minimum as follows:

- recommended drill speeds:
2.0 mm Ø: approx. 800 rpm, 2.75 mm Ø: approx. 700 rpm,
3.0 mm Ø and 3.4 mm Ø: approx. 500 rpm, 3.7 mm Ø,
3.9 mm Ø, 4.25 mm Ø, 4.9 mm Ø and 5.25 mm Ø:
approx. 400 rpm
- intermittent drilling and cutting using sharp drills and cutters (use 10 times only)
- cool the instrument externally and internally with sterile physiological saline (NaCl) or Ringer's solution.

Important

To determine the optimum implant length and prevent nerve injuries, check in advance that adequate bone is available (X-ray).

Placing the implant after preparing the implant bone area

We advise to study our detailed surgery documentation CD/System overview thoroughly.

Removing the insertion post

In order to remove the insertion post (red, blue, brown), pull off the plastic implant carrier and loosen the vertical screw (blue) with the hexagonal screwdriver 1.2 mm. Dispose of the vertical screw. In cases with concomitant augmentation or where the primary stability is low, secure the implant against rotation with help of a counter wrench. Afterwards, the insertion aid can be lifted from the implant easily.

Closing the soft tissue

The cover/healing aid must be screwed onto the implant before the wound is closed with atraumatic suturing material.

Recovering the implant

Once the individual healing period (minimum of three months for the mandible and six months for the maxilla) has elapsed, the implant should be recovered and the cover/healing screw removed. In case of augmentation techniques, the healing period must be prolonged accordingly. The interior of the implant is then cleaned and the correct height of gingiva former for the oral mucosa is screwed in. Make sure that the gingiva former fits properly! The restoration should not be fitted to the implant until the soft tissue has healed fully and shows no irritations.

Impression taking

Once the gingiva former has been removed, the impression post should be screwed into the implant. When taking the impression, the impression material is first syringed around the impression post. Once the impression material has cured, loosen the screw, remove the impression from the mouth together with the impression post and close the implant again with the gingiva former. Carefully screw an implant analog into the impression post and cast the model – preferably including a gingival mask. Once the cast has set, remove the impression post, respectively transfer coping, and replace it with the required abutment.

Laboratory components

When selecting the abutments, the angulation of the implant and the thickness of the gingiva must be taken into account. Those areas of the abutment which are in contact with the implant may not be blasted with abrasives or treated in any other manner.

Documentation

Every IMPLA implant can be traced via its lot number. The label on the outer blister package includes all the data for recording in the patient's file. Any claims can only be asserted when you provide us with the requested information, indicating the above mentioned data.

Packaging and sterility

IMPLA implants are packed sterile in a special glass container. The chosen implant is removed from the sterile container with the insertion aid only directly before insertion and is inserted into the prepared bone cavity without touching the implant.

A special attachment on the bottom of the implant container holds the sealing, covering and healing screw which can be easily removed with the screwdriver 1.2 mm wide hexagon wrench key.

The outer cardboard box shows the user important data such as implant length, implant diameter, article no., lot and CE certification.

General information

- All implants are for single use only.
- Each is packed sterilely in a closed container, please note the expiration date of sterility.
- The implant is sterile (beta sterilization) as long as the packing remains unopened and undamaged. A subsequent sterilization by the user is not permitted by the manufacturer. A chemical sterilization of the implant (chemiclave) may not be carried out under any circumstances, as this could damage the implant surface. In case of damaged packaging, the implant is no longer sterile and may not be used anymore.
- Do not touch the implant.
- Store the sealed and undamaged implant container in a dry place at room temperature. Avoid direct light incidence.
- System errors due to mixing up instruments and implants should be avoided. Please always check the color codings and labelings.

- Thermal damages of the jaw bone by rotating instruments (such as pilot drills and extension drills) must be avoided.

- Please find information on the preoperative diagnostics, indications and contraindications, risk factors, post-operative treatment of the implant as well as patient information sheets, the „Informed consent“ form and surgery record form „surgery log file“ in our IMPLA implants instruction manual.

Torque of screw connections

10 Ncm horizontal screw connections
20 Ncm vertical screw connections
30 Ncm insertion of the implant, Locator Abutment screw connections

CE 0297

SCHÜTZ  **DENTAL GROUP**

Schütz Dental GmbH
Dieselstr. 5-6 · 61191 Rosbach/Germany
Telefon: +49 (0) 6003 814-0 · Telefax: +49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de