Alphalink

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Directions for use

Alphalink is used to adhere metal components directly on the working model, ensuring an absolutely passive fit.

A high proportion of inorganic fillers guarantees a high final hardness producing the stability required for laboratory preparation.

Indications

- Fixing telescopic or conical crowns to chrome cobalt frameworks.
- Fixing pre-fabricated attachments to chrome cobalt frameworks.

Constructing the joint

The areas to be adhered must be designed so as to be retentive and strong.

When the secondary component is placed on the primary component, an even, thin gap must be present.

This gap, which ensures that a passive fit is obtained, is then filled with Alphalink.

When adhering pre-fabricated attachments, undercuts must be created in the joint areas and the manufacturer's instructions heeded.

Preparing the joint areas

After ensuring that both the primary and secondary components fit passively, create undercuts in them and sandblast with aluminium oxide (250 µm).

After steam cleaning, condition the metal surfaces with Sebond Smart metal bonding system.

Application

Alphalink is available in the user-friendly dual-compartment cartridge as well as in two individual syringes.

Please proceed as follows:

- Place a disposable mixing tip onto the **dual-compartment cartridge**. Discard the first bit of mixed material discharged from the tip.
- Press Alphalink from the automix tip and apply the desired amount directly onto the restoration. Spread evenly.
- As the material will cure inside the used mixing tip, leave the tip on the cartridge as a closure until the next use.

- By using the automix cartridge with the mixing tip, Alphalink is always mixed in the ideal ratio of 1 : 1.

When using **individual syringes**, please proceed as follows:

Two equal lengths of components A and B are dispensed from both syringes by pressing the plunger and are manually mixed on a mixing slab to a homogeneous consistency using a spatula.

Leave the mixing cannula as a protective cap on the needle until next use.

Working time after starting mixing:

Approx. 3 min

Setting commences after approx. 4 1/2 min

Setting completed after approx. 8 min

These times assume a room temperature of 22 °C. Higher temperatures accelerate the setting reaction and lower temperatures retard it.

Application

Alphalink applied to matrix and patrix. The metal components should be repositioned on the model and fixed into place.

Ensure that the areas to be bonded are completely filled with Alphalink and that there are not any air bubbles or contamination.

Important

The sections being joined must not be moved until the Alphalink has set completely.

Please note

- Alphalink is dimensionally stable up to 120 °C.
- Alphalink bonds to commercially available denture materials.
- Alphalink is easy to grind.

Storage

If refrigerated at 3-9 °C. Do not use after the expiry date.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies)

or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications/interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervisi-

on. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient’s mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Troubleshooting

Working time too short	– High ambient temperatures (e. g. summer heat or central heating) decrease the working time.
Working time too long	– Alf the material is used directly after taking it from the refridgerator, curing will be delayed.
Material does not cure sufficiently	– Substances containing eugenol or oil of wintergreen prevent curing. – Material components are insufficiently mixed; only use original mixing tip.
Material is cured	– After use, the cap was replaced instead of leaving the mixing tip on the cartridge (contamination of component A +B).
Construction does not fit	– Construction was moved before the material had cured. – Material is already too firm: Please join fewer parts at a time.

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