FuturaBasic Hot

CE 0297

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Tel.: +49 (0) 6003 814-0 · Fax: +49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de

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Instructions for use

Application:

Fabrication of partial and full dentures using the flasking and pressing technique

Characteristics:

- Colour stable
- Easy to use

Mixing ratio:

2.5 parts powder by weight

1 part liquid by weight

Approx. **25 g powder** and **10 g liquid** are required for a full denture.

Technique:

- Coat the plaster surfaces with FuturaSep Plus alginate separating agent (Schütz Dental).
- After adding the powder to the liquid, spatulate FuturaBasic Hot thoroughly to a homogeneous mixture in a mixing vessel for 30 secs.
- Allow a doughing time of approx. 15 minutes. The acrylic dough then has a working time of approx. 30 minutes.
- At the end of the doughing time knead the acrylic dough thoroughly for approx. 30 seconds.
- When the FuturaBasic Hot becomes viscoelastic, pack it into the prepared flask and press.
- Hold the pressure for approx. 15 minutes.
- Secure the flask in a clamp for polymerisation.
- The denture teeth should be roughened, provided with undercuts and coated with monomer or HS Cross Liquid (Schütz Dental) to ensure they are securely retained.

Polymerisation:

Polymerise FuturaBasic Hot in the closed flask with the clamp for 20 minutes in boiling water. Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

Repairs:

FuturaBasic Cold or FuturaSelf can be used for repairs. Before applying the acrylic dough, the fractured surfaces should be roughened, cleaned and coated with monomer. Polymerise in a pressure vessel for 15 minutes at a pressure of 2-5 bar with a water temperature of 45° C.

Storage:

Store at 10-25° C. Do not use after the expiry date.

Important:

Should the room temperature deviate from 23° C: higher temperatures reduce the working time and lower temperatures extend the working time.

Presentation:

Powder in 1.000 g, 3.000 g and 12 kg

Liquid in 1.000 ml

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare.

Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-Indications/interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Caution:

Avoid skin contact with unpolymerised material and inhalation of monomer fumes.

Disposal:

Powder can be disposed of as domestic waste. Liquid should be disposed of as hazardous waste.

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