# Tizian Blank Zirkon-verstärktes Komposit

---

CE 0297

Schütz Dental GmbH · Dieselstr. 5 - 6 · 61191 Rosbach/Germany  
Tel.: +49 (0) 6003 814-0 · Fax: +49 (0) 6003 814-906  
[info@schuetz-dental.de](mailto:info@schuetz-dental.de) · [www.schuetz-dental.de](http://www.schuetz-dental.de)

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## Instructions

### Description

Tizian Blanks Zirconia Reinforced Composite are made of a composite material reinforced with zirconium (PMMA-based, 20 % anorganic filler). The blanks can be processed with CAD/CAM technology. Optionally, they can be veneered with a composite or with a monolithic material. Tizian Blanks Zirconia Reinforced Composite are available in several different tooth shades.

### Indications:

#### 1. Long-term temporaries (for up to 2 years)

- a) Up to 16-unit anterior and posterior crowns and bridges, fully anatomic or reduced.
- b) Position no more than 2 units in the posterior area.
- c) If the bridge is a free-end construction, add max. one posterior attachment.

#### 2. Smaller permanent restorations\*

- a) individual crowns
- b) inlays
- c) onlays
- d) veneers

#### e) bridges with a maximum of 3 units

\*It is highly recommended to veneer such restorations with an incisal material.

If the restoration is not veneered with a composite material, a high lustre polish (smooth surface) is indispensable.

### Material properties

- Composition: > 78 % PMMA, 20 % inorganic filler, < 2 % anorganic pigments
- E-module: 3050 MPa
- Flexural strength: 72 MPa

- Hardness: 196 MPa
- Water solubility: 1.6 µg/mm<sup>3</sup>
- Residual monomer: < 1.0 %

### Framework design

- Minimum wall thickness:
  - 0.6 mm (cervical)
  - 1.0 mm (occlusal)
  - > 0.2 mm (crown margin)
- Diameter of bridge units: max. 10 mm
- Profile of connectors:
  - > 10 mm<sup>2</sup> (anterior)
  - > 15 mm<sup>2</sup> (posterior)

**Important:** To increase the stability of the construction, choose the greatest possible height for the connectors.

### Milling procedures

- The optimal ambient temperature for milling this material is 20 °C – 25 °C (68 °F – 77 °F).
- Please adhere to the instructions of the milling system. We recommend to use a **solid 5-axes milling system** to ensure an excellent fit of the restorations.
- Tizian Blanks Zirconia Reinforced Composite are suited to dry milling.
- All milling **tools** must be suited to processing composites. We recommend coated radius cutters which are conventionally used for composites.
- The choice of **milling parameters** strongly depends on the milling system, on the CAM-software and on the milling geometries used. If you are not certain which parameters to use, please use your parameters for composites.

### Detaching the restoration

Carefully detach the finished restoration from the blank with a jig cutter.

## **Finishing / polishing**

- Use clean cross-cut hardmetal mills to finish the restorations.
- Use a low number of revolutions and a minimum of pressure to avoid thermal damage. Ensure sufficient cooling.
- Polish with conventional pastes for polishing composites (e. g. dialog Vario polish), goat bristle brushes and cotton polishers.

### **Important:**

- Do not overheat the material when using mills, steam cleaners or similar.
- The wall thickness of the frameworks must not go below a limit of 0.4 mm.

## **Cleaning the frameworks**

- Alcohol in high concentrations will generally cause cracks on PMMA, the so-called crazing effect. Do not use cleaners or disinfectants which contain alcohol.
- Do not overheat the surface of the restoration while steam-cleaning it.
- Alternative possible cleaning methods are mechanical procedures, using an ultrasonic unit or alcohol-free disinfectants with quaternary ammonium compounds.

## **Veneering**

After suitable surface conditioning, Tizian Blanks Zirconia Reinforced Composite can be veneered with all conventional veneering composites:

### **Recommended procedure**

1. Condition the surfaces by sandblasting with aluminium oxide (max. 110 µm and 2 bar, obtuse angle). Ensure that the surfaces are clean, dry and grease-free, before you continue processing the material.
2. Apply a suitable bonding agent, e. g. dialog Bonding Fluid (Schütz Dental).

3. Veneer conventionally, e. g. with dialog Vario Composite Paste and/or Flow. Observe the instructions of the composite's manufacturer.

**Important:** When preparing a vestibular veneer, it is necessary to apply a protection on the masticatory edge.

### **Tip:**

If you are using the Schütz Dental veneering system, dialog Occlusal masses can be applied directly to the framework, without a foundation of paste or flow material. To obtain an even bond in such cases, we recommend the use of Sebond Grip.

When using a different veneering system, the incisal mass should not be applied directly to the framework. This could lead to chipping, depending on the material used.

## **Fixation**

The material must only be fixed by adhesion with a conventional temporary or permanent fixation composite. Please adhere to the instructions of the fixation composite's manufacturer.

### **Note**

Do not use cements of any kind.

### **Safety note**

During processing of „Tizian Blanks Zirconia Reinforced Composite“, dusts are released which can damage the lungs, and irritate the skin and the eyes. Therefore, please only process the material while running an adequate extractor system.

Please adhere to the following safety rules:

- Process the blanks only underneath a local extractor system.
- Ensure that product dust will not come in touch with the eyes or the mucous membranes.

- Do not inhale the dust.
- Wear a face mask and protective goggles.

### **Storage**

Store in a dry place at temperatures between 10 °C and 25 °C (50 °F and 77 °F). Keep away from direct sunlight and moisture.

### **Disposal**

The material can be disposed of with normal household waste.

### **Guarantee**

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines.

### **Shelf-life**

The maximum shelf-life is printed on the label. Do not use after the shelf-life has expired.

### **Side effects**

Undesirable side effects of this medical device are to be expected only in extremely rare cases as long as the product is handled and used properly. However, immune reactions (e.g. allergies) or local discomfort cannot be excluded completely. If you become aware of any side effects, please notify us, even in case of doubt.

### **Contra-indications / interactions**

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

### **Note**

Please ensure that all above information is readily available to the attending dentist in case this product is used. Please also note the material safety data sheets available.

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