AlphaForm

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Directions for use

AlphaForm – For duplication of master models to process Chrome Cobalt and acrylics.

AlphaForm – Especially suitable for production of study and teaching models.

AlphaForm is very precise (0,06% contraction)
– Compatible with all investment, plaster and acrylics.

AlphaForm – has a long shelf life.

AlphaForm – Multiple models can be made from one mould.

AlphaForm – Simple, quick and safe to process.

AlphaForm – has resistance to solvents such as alcohol etc.

AlphaForm – has good flowing properties and penetrate every area of a complicated model.

Instructions for use

1. The master model must be dry (not water soaked) and placed in a normal duplicating flask.
2. The required amount of AlphaForm base and hardener are mixed in the ratio 1:1 for approx. 1–2 minutes until a homogenous consistency is achieved.

To ensure optimum use, we recommend our automatic mixing and dosing device „Alphamat“.

3. Pour the mixture into the duplicating flask.
4. Working time 4–5 minutes.
5. Deflask after 45 minutes.

Setting time

Within 45 minutes at normal room temperature. These times will be altered by changes in room temperature: i.e. higher temperature quicker setting.

Note

Bubbles in the surface of the mould will be avoided if the AlphaForm is poured slowly into the mould, and if the mould is first sprayed with a silicone releasing agent from Schütz Dental.

It is important that the master model is perfectly dry.

Fine parts of the model must be brushed or sprayed.

Fault Analysis

Manual mixing:

Fault	Cause	Solution
Material does not solidify	insufficient stirring or no stirring	observe mixing ratio 1:1 stir carefully
Material solidifies unevenly	insufficient stirring	stir carefully

Alphamat:

Fault	Cause	Solution
Material does not solidify / Material solidifies unevenly	<ul style="list-style-type: none"> • no / wrong mixing cannula • insufficient stirring 	<ul style="list-style-type: none"> • use correct mixing cannula • take care that there is an even output of both components • check level of content in reservoir and refill if necessary

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