# FuturaPress LT

**MANI**® **SCHÜTZ**  
**DENTAL**

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany  
Telefon +49 (0) 6003 814-0 · Telefax +49 (0) 6003 814-900  
info@schuetz-dental.de · www.schuetz-dental.de

**CE 0297**

Mandler 04/2017 - 1.000 - #500003

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## Instructions

### Indication:

For completion of total dentures, partial dentures and frameworks, for repairs, extensions and relinings in the flow technique.

### Properties:

- colorfast
- high mechanical values
- easy to process
- long processing time

### Mixing:

10 parts by weight of powder  
7 parts by weight of fluid

### Processing:

- Insulate the plaster model with the alginate insulation FuturaSep Plus (Schütz Dental).
- After adding the plastic powder to the fluid and letting it thicken for about 15 sec., stir FuturaPress LT to a homogeneous mixture in a beaker.
- During the pouring phase (2-3 min.), the FuturaPress LT paste can be filled into a silicone index, for example. After the plastic phase starts, it will no longer flow out of the index and can be shaped.
- To prevent plastic teeth from breaking off, they must be roughened, as usual for PMMA plastics, have undercuts made, and be moistened with monomer or Cross Liquid (Schütz Dental).
- After max. 8 min. (at an ambient temperature of 23 °C), FuturaPress LT must be put into a pressure polymerization device (e. g. Futura-mat from Schütz Dental) for polymerization.

### Important:

In case of deviations from a 23 °C ambient temperature – higher temperatures shorten, lower temperatures lengthen the processing times. Do not use after the expiration date.

### Polymerization:

The polymerization time in a pressure polymerization device is 15 min. at a pressure of 2 - 4 bar and water temperature of 45 °C.

### Note on repairs:

Repairs can be made with FuturaPress LT, FuturaPress N, and FuturaSelf for example. Before applying the prepared plastic paste, breakage areas must be roughened, cleaned,

and moistened with monomer. The polymerization time in a pressure polymerization device is 15 min. at a pressure of 2 - 4 bar and water temperature of 45 °C.

### Storage:

Keep at 10 - 25 °C.

### Available as:

Powder in 500 g, 1000 g, 12 kg.  
Fluid in 500 ml, 1000 ml.

### Hazard note:

Avoid prolonged skin contact with unpolymerized material and inhalation of monomer vapours.

Once the container has been opened: Do not use after the expiry date.

### Disposal:

The powder can be disposed of with household waste. Dispose of the liquid in a special refuse tip.

### Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

### Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and cross-reactions of the product with other materials already in the patient's mouth before using the product.

### Hazard warning:

Avoid skin contact with unpolymerized material for a longer period of time and inhalation of monomer vapours.

Open containers: Do not use after the expiration date.

## Troubleshooting

Problem	Cause	Remedy
Acrylic does not fully cure	Incorrect mixing ratio	<ul style="list-style-type: none"> <li>– Adhere to the instructions for use</li> <li>– Adhere to the mixing ratio</li> </ul>
Acrylic is whitish on the fitting surface	<ul style="list-style-type: none"> <li>– incorrect or poor application of separating agent</li> <li>– model not soaked</li> </ul>	<ul style="list-style-type: none"> <li>– Apply ample separating agent</li> <li>– Check the separating agent if necessary</li> <li>– Soak the model for approx. 10 min. in lukewarm water</li> </ul>
Acrylic forms bubbles on the fitting surface	Model not soaked	<ul style="list-style-type: none"> <li>– Soak the model for approx. 10 min. in lukewarm water</li> </ul>
Acrylic surface white / porous	Recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Color differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic / denture has whitish streaks	Incorrectly mixed / too much powder used	Thoroughly spatulate the acrylic dough / adhere to the mixing ratio
Acrylic sets too fast	Temperatures too high Working stages too long	Observe temperatures / working times Storage: 10 - 25 °C

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