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### Ceramic-Composite Repair System

#### Indication

The combination of the ReVeneer Primer and Base represents a fast and universally applicable system for intraoral and extraoral repair of fractured ceramic blending without acid cauterization.

#### Preparation of Surface Fractures

Smooth the damaged surface and roughen it with a fine grit diamond (< 50 micron) or intraoral air polisher before cleaning off any adhering material with oil-free compressed air.

#### Please note

The use of an intraoral air polisher ( $\text{Al}_2\text{O}_3 / < 120$  micron) further enhances the bond strength.

#### ReVeneer Primer

##### Application

In combination with ReVeneer Base, ReVeneer Primer serves as an adhesive agent for ceramic/metal and composites.

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#### Troubleshooting-Liste

## Troubleshooting List

Errors	Cause	Remedy
Blending does not hold after repair	<ul style="list-style-type: none"><li>- Surface not sufficiently roughened, cleaned and dried</li><li>- ReVeneer Opaque or Base not entirely hardened</li><li>- Faulty frame dimensions</li><li>- The repaired composite was not placed out of occlusion.</li></ul>	<ul style="list-style-type: none"><li>- rough out frame surfaces and remove all grease layers, moisture and milling dust</li><li>- see following points</li><li>- new dimensioning of frame</li><li>- In areas exposed to severe masticatory loading, trim the composite to prevent occlusal contact.</li></ul>
ReVeneer Opaque not completely hardened	<ul style="list-style-type: none"><li>- Opaque applied too thickly</li><li>- Single coats not polymerized separately</li><li>- Insufficient light output by polymerization lamp</li></ul>	<ul style="list-style-type: none"><li>- thinner opaque application</li><li>- separate polymerization of individual coats</li><li>- check light output and eventual change of light source</li></ul>
ReVeneer Base not completely hardened (paste still green after light hardening)	<ul style="list-style-type: none"><li>- Base applied too thickly</li><li>- Insufficient light output by polymerization lamp</li></ul>	<ul style="list-style-type: none"><li>- apply thinner coat of ReVeneer Base</li><li>- check light output and eventual change of light source</li></ul>
ReVeneer Base in the needle appears to be too hard and solid	<ul style="list-style-type: none"><li>- Material not correctly sealed, paste pre-polymerized</li><li>- Material not warmed to room temperature after removal from the refrigerator</li></ul>	<ul style="list-style-type: none"><li>- replace cap on needle after each material withdrawal</li><li>- allow ReVeneer Base to warm to room temperature before use</li></ul>
ReVeneer Base in the needle appears too fluid		<ul style="list-style-type: none"><li>- observe the storage temperature</li></ul>

#### Warranty

Our technical user tips, regardless of whether they are provided verbally, in writing or through practical guidance, stem from our own experiences and can therefore be regarded only as guidelines.

Our products undergo continuous development.

Therefore, we reserve the right to possible future changes.

#### Indications

Please observe the existing relevant safety datasheets while processing.

#### Side Effects

With appropriate processing and use of this medical product, undesired side effects are rarely experienced. Immune reactions (such as allergies) or local irritations cannot be ruled out entirely in principle. Please inform us of any undesired side effects you may have experienced – even when in doubt.

#### Adverse Signs / Interactions

Do not use or use only under strict supervision of the attending physician / dentist if patient suffers from hypersensitivity. In this instance, the composition of our medicinal product is available upon request.

Known cross-reactions or interactions of this medicinal product with other materials already in the mouth must be taken into account by the dentist prior to use.

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Fotopolimerizzatore	ReVeneer Opacizzante	ReVeneer Base
Spektra 2000 (Schütz Dental)	180 sec.	180 sec.
Spektra LED (Schütz Dental)	50 sec.	20 sec.
Uni-XS (Heraeus Kulzer)	180 sec.	180 sec.
Optilux 501 (Kerr)	60 sec.	40 sec.
CapoLux (Schütz Dental)	60 sec.	40 sec.
Translux CL (Heraeus Kulzer)	60 sec.	40 sec.

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