Flexor CC

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Instructions for use

The Flexor CC Combi pack is for soft lining dentures to make them much more comfortable to wear. Flexor CC is an automixing, permanently soft direct relining material, based on addition-curing silicone and supplied in cartridges (AUTOMIX).

This system comprises Flexor CC (silicone), Flexor CC primer and Flexor CC glaze. Flexor CC primer bonds Flexor CC to the denture and Flexor CC glaze (base and catalyst) smoothes areas of soft liner roughened by trimming.

Flexor CC contains modified polydimethyl siloxane with Pt catalyst and the primer contains dichloromethane (which evaporates and does not contact the patient). In case of allergies against these substances, do not use this product.

Applications and properties

- Suitable for all PMMA acrylics
- Permanently soft relines are quickly and easily fabricated directly in the patient's mouth or indirectly in the laboratory
- Avoids sore spots on the upper and lower jaws
- Seals the palatal vibrating line in cases where the denture is difficult to retain or the alveolar ridge is flabby
- Creating an elastic periphery on upper dentures to increase the adhesion and enhance retention
- Covering sharp alveolar ridges with an elastic material
- In implantology for sealing implants and overdentures during the healing period

Directions

1. Preparing the denture

The existing relining material should usually be removed completely. Trim the denture base to reduce it by 1-2 mm. Trim a 1 mm ledge into the vestibular and lingual/dorsal peripheries at a height of 2-3 mm. Round off the denture periphery. Roughen the areas to receive the relining material.

Once the denture has been prepared, clean it thoroughly and use alcohol (ethanol or isopropyl alcohol) to remove all traces of fat.

2. Applying Flexor CC primer to bond the material

Apply a uniform coat of Flexor CC primer to the cleaned and trimmed surfaces to which the soft lining material is to be bonded, ensuring that the entire surface is coated (beyond the peripheries). Allow the primer to dry for approximately 1 minute. Close the bottle immediately after use. Saliva, moisture or monomer liquid prevent Flexor CC bonding to the denture.

3. Preparing the dispensing gun and cartridge

Fit the double cartridge (AUTOMIX) to the dis-

persing gun. Remove the cap from the AUTOMIX double cartridge and dispose of it. **The cap must not be used again.** Extrude the material gently until it flows out of both cartridge outlets uniformly. Fit the mixing tip to the cartridge and extrude the required amount of material. After use, leave the mixing tip on the cartridge to seal it.

4. Relining the denture

4.1 Relining a denture directly in the patient's mouth

Flexor CC (AUTOMIX) is mixed automatically and fully with a dispensing gun. Once the material has been mixed, apply it uniformly to the vestibular and lingual/dorsal peripheries. Then coat (at least 2 mm thick) the remaining area of the denture base. **It is important to coat all areas of the denture in contact with the soft tissues.**

The coat of soft lining material must be applied **within 1 minute** (working time). Then place the denture in the patient's mouth. The patient should then close his mouth, applying normal occlusal loading, and carry out all functional movements. The material sets in **5-6 minutes**. The soft lined denture is then removed from the mouth, trimmed and glazed (refer to sections 5 and 6).

4.2 Indirect method in the dental laboratory

Flask the denture as usual, with an opposing model. Remove the denture from the flask and prepare it as described in sections 1 and 2. Coat the model with a plaster/acrylic separating agent (e.g. Futura Sep Plus). Use the double cartridge and dispensing gun to apply a uniform coat of Flexor CC to the surface of the denture to receive the soft lining. This coat should be approximately 1 mm thicker than the acrylic trimmed off previously. Close the flask with the opposing model, place the flask in a clamp and insert it into a pressure vessel for 10 minutes at 40-45 °C. Then deflask the denture and trim/glaze it as described in sections 5 and 6.

Please note: These times apply at 23 °C. Lower temperatures (e.g. after refrigerating the material) prolong, and higher temperatures curtail, the times. Do not use together with condensation-curing silicone.

5. Trimming the denture

Remove the bulk of the excess material with sharp scissors or a scalpel. Then use a silicone trimmer, e.g. the enclosed polishing wheel, and exert gentle pressure to smooth the edges of the cuts.

6. Glazing the denture

The enclosed Flexor CC glaze is used to smooth the functional peripheries after trimming. 5 drops each of Flexor CC glaze base and Flexor CC glaze

catalyst are placed in a mixing well or on a mixing pad, mixed with a spatula and brushed on. Allow the glazed surfaces to dry for 5 minutes at room temperature.

7. Recommended cleaning procedure for patients
Scrub and rinse the relined denture thoroughly under running water several times per day using a toothbrush and neutral soap. Cleaning solution (e. g. in a Microclean): up to 10 minutes.

8. Technical data – Flexor CC

Working time: Approx. 1 minute
Setting time: Approx. 5:30 minutes
(incl. working time)
Recommend intraoral setting time:
Approx. 5 - 6 minutes
Polymerization shrinkage: Less than 0.1%
Permanent set: Less than 1%
Shore a hardness: Approx. 33

9. Flexor CC products

Flexor CC set (50 ml cartridge, 4 ml primer, 2 x 5 ml glaze, 10 mixing tips, 10 brushes, 2 mixing wells, 1 polishing wheel)
Order No.: 250030

Flexor CC cartridge 50 ml
Order No.: 250031

Flexor CC glaze 2 x 5 ml
(base and catalyst)
Order No.: 250033

Mixing tips, 50 pieces
Order No.: 212910

Polishing wheels, 25 pieces
Order No.: 250024

Dispensing gun
Order No.: 212914

10. Storage

Store in a dry place at 10 - 25 °C.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.