

# Capo Bulk Fill

Quick, Easy and Safe



uberVo

DENTAL B.V.

Importeur | distributeur  
Schütz Dental NL

# Capo Bulk Fill – Quick, Easy and Safe

Work quicker and easier than ever before without compromising the physical properties.

Saves time with layers of up to 4 mm.

Capo Bulk Fill is a light-curing composite for posterior restorations in the direct filling technique and for the “Bulk Fill”-technique. This quick and easy method is suitable for layers with a thickness up to 4 mm. The stable flow material adapts really well to the bottom and walls of the cavity.

## Indications

- Fillings with layer thicknesses up to 4 mm in Black's class I cavities.
- Base fillings/liners with layer thicknesses up to 4 mm in Black's class II cavities. For the final layer in a class II cavity only use composites which are suited to filling Black's class II cavities.
- Cavity lining – as a first (bottom) layer in Black's class I and II cavities.

## Benefits

- No additional layers are necessary to cover class I cavities. Saves time.
- Quick and easy to work with, suitable for layers with a thickness of up to 4 mm
- Minimal shrinkage provides for a very safe application
- Outstanding physical properties (filler content approx. 77%)
- High radiopacity (210% Al) ensures a dependable diagnosis between the filling material and the natural tooth substance
- Compatible with most bonding and composite systems
- Biocompatible: 100% free of Bis-GMA, HEMA and TEG-DMA
- Universal dentin shade provides a natural look
- Excellent flow behavior
- Easy to polish

## Physical data

	Mean value
3-point-strength (MPa)	134
Young's modulus (MPa)	11.090
Vickers hardness (MPa)	737
Compressive strength (MPa)	381
Radiopacity	210 % Al
Filler content	77 % by weight



art. no. 232060

2 g syringe including 10 application tips.

Universal shade adapts to the natural tooth color.


## Expert opinion Dr. Schulze



Dr. med. dent. Ellen Schulze,  
Wimmelburg/Germany


*“Capo Bulk Fill is very easy and quick to shape and to polish. Its great advantage is its flow behavior in undercuts and in small approximal preparations. The product offers a great cost/performance ratio. I am always happy to recommend the material to my colleagues.”*

## Schütz Dental



Capo Bulk Fill (1) "one step filling". \*

## Competitors



Conventional bulk fill material (2) with occlusal cover (3).

\* Applies to the cavity classes I, III, IV and V.

## Capo Bulk Fill Compules

Benefit from both the outstanding physical properties of the Bulk Fill composite and the convincing advantages of the compules.

Work quicker, safer and more efficiently.

Advantages of the compules:


- High output quantity
- Hardly any loss of material
- First-rate surface quality
- Especially hygienic use
- Cap is easy to remove



# Physical properties

## Scientifically tested by renowned institutions


### Physical tests



Source: Fraunhofer Institut für Keramische Technologien und Systeme, Dresden/Germany  
September 2012

### Volume shrinkage test


Measuring the volume shrinkage at 37°C by initializing with blue-light.



Source: InnoMat GmbH, Teltow/Germany  
07.09.2012

### Cytotoxicity test

"In this study under the given conditions no leachable substances were released in cytotoxic concentrations from the test items."



The test showed that Capo Bulk Fill has no cytotoxic effects.

Source: BSL BIOSERVICE Scientific Laboratories GmbH Planegg/Germany  
Study No.: 123515  
10.07.2012

# The perfect fit

## Capo Natural

Ultra-fine micro hybrid composite with improved abrasion characteristics. Available in different tooth shades. High biocompatibility.



art. no. 232055

## Capo Bond Set

Multifunctional, light curing 1-component primer and bonding system, ethanolbased. Excellent bonding strength on dentin and enamel.



art. no. 234014

## Capo Universal

Has the same great features as Capo Natural. Available in one universal color for posterior restorations. High brilliance.



art. no. 232054

## Capo Etch Set

Safe and fast etching of dentin and enamel. High viscosity and excellent contrast to natural teeth, which makes the product easy and safe to apply.



art. no. 234013

## Capo Slow Flow

Combination of the easy handling of flowable composite with the exceptional physical properties of paste. Free of Bis-GMA. Available in 8 colors.



e.g. art. no. 232062

## Capo Slow Flow Compules

The easy handling of a flowable composite and the outstanding physical properties of a paste – now available in compules!



art. no. 232112-232115

## ReVeneer Adhesive Set

Ceramic repair system consisting of primer, repair paste and composite. Aesthetic, lasting, safe. Free of acid and silane.



art. no. 232205



**Contactgegevens:**

uberVo dental b.v.  
Electronweg 10d  
1627 LB Hoorn NH NL  
T. +31(0)229 213736  
info@ubervo.nl  
www.ubervo.nl

importeur | distributeur

**SCHÜTZ**  **DENTAL**  
Micerium Group

**uberVo**  
DENTAL B.V.